

# DM-120 Digital Micro Centrifuge

Laboratory Bench-top Centrifuge

Centrifuge



## Features

- Auto-balancing system to extend centrifuge life as well as to reduce noise and vibration
- Short run option with the press of a button
- Precise digital timer with 30 sec increments
- Safety interlock to prevent lid from being opened during operation
- Standard rotor for micro centrifuge tubes, optional hematocrit rotor available

## Specifications

Item	Specification
Speed	500—12,000 RPM
Timer	0 - 15 min digital timer w/ 30 sec increments
Display	LED digital displays rpm and remaining time
Power	110V 60Hz or 220V 50Hz
Weight	32 lb
Size (L x W x H)	290 x 400 x 260 mm
Construct	Metal Housing
Micro centrifuge tube Rotor <i>(Standard)</i>	45° angled rotor 24 x (1.5~2.0) ml 500-12,000 RPM Max. Rcf 16,096 xg
Hematocrit Rotor <i>(Optional)</i>	Disc Rotor 24 x capillary tubes (dam 1.5mm x 75mm) 500-12,000 RPM Max. Rcf 14,648 xg Standard aluminum percentage reader with purchase of hematocrit rotor



Micro centrifuge tube rotor



Hematocrit Rotor



Aluminum percentage reader